

FUSION BOND 7

Light Curing Self Etching

One Component Bonding Adhesive

PREVESTDenPro®
THE FUTURE OF DENTISTRY

Product Description : Fusion Bond 7 combines the ingredients needed for etching, priming, and bonding into a single adhesive, thereby eliminating separate etching and priming steps of the bonding process. Fusion Bond 7 adhesive can be used for the bonding of both direct and indirect restorations.

Composition : Hema Phosphate, Urethane Dimethacrylate, Triethylene Glycol Dimethacrylate, Bisphenol A Diglycidyl Methacrylate, 2-Hydroxy Ethyl Methacrylate, Water, Tertiary Butanol, Photoinitiators, Stabilizer.

Indications :

Direct Applications :

1. Light-cured composite and compomer restorations.
2. Composite/ceramic/ metal repairs.
3. Cavity sealing for amalgam restorations.
4. Light-cured or dual-cured core build-ups (always light-cure an incremental build-up of core material).

Indirect Applications :

1. Veneers
2. Porcelain, composite, and metal-based inlays, onlays, crowns, bridges.
3. Endodontic posts.
4. Cavity sealing as a pretreatment for indirect restorations.

Instruction for use : Bonding Composites/Compomers to Dentin/Enamel Veneers, crowns, bridges, inlays, onlays and metal-based restorations (used in conjunction with a luting cement)

1. Prepare cavity. Wash thoroughly with water spray and air dry (Do not desiccate). **Note :** For optimal adhesion, always prep (bevel) all enamel and/or margins prior to application of Fusion Bond 7 adhesive.

2. Shake Fusion Bond 7 adhesive bottle briefly. **Note :** Upon first use, shake Fusion Bond 7 adhesive bottle vigorously for 20 - 30 seconds.

3. Dispense 2-3 drops of Fusion Bond 7 adhesive into a clean well. Replace cap immediately after dispensing.

4. Using the disposable applicator brush, apply a generous amount of Fusion Bond 7 adhesive to the enamel/dentin surface. Scrub the surface with a brushing motion for 20 seconds.

5. Apply a second application of Fusion Bond 7 adhesive with a brushing motion for 20 seconds.

6. Dry the adhesive with gentle air first and then medium air for at least 5 seconds with oil-free air.

7. Light cure for 20-30 seconds.

8. Apply composite according to manufacturer's instructions for use.

Instructions for bonding Indirect Restorations :

1. Clean the preparation. Wash thoroughly with water spray and air dry (Do not desiccate).

2. Shake Fusion Bond 7 adhesive bottle briefly. **Note :** Upon first use, shake Fusion Bond 7 adhesive bottle vigorously for 10 sec.

3. Bottle Delivery : Dispense 2-3 drops of Fusion Bond 7 adhesive into a clean well. Replace cap immediately after dispensing.

4. Using the disposable applicator brush, apply a generous amount of Fusion Bond 7 adhesive to the enamel/dentin surface. Scrub the surface with a brushing motion for 20 seconds.

5. Apply a second application of Fusion Bond 7 adhesive with a brushing motion for 20 seconds.

6. Dry the adhesive with gentle air first and then medium air for at least 5 seconds with oil-free air. Avoid pooling of adhesive before light curing.

7. Light cure for 20-30 seconds.

Contraindications :

1. Do not use Fusion Bond 7 adhesive if a patient is known to be allergic to any of the material's ingredients or if the stipulated working technique cannot be used.

2. Fusion Bond 7 adhesive must not be used in conjunction with purely self curing composite materials. Dual-curing materials must be light activated.

3. Do not use Fusion Bond 7 adhesive for direct pulp capping.

Warning : Fusion Bond 7 adhesive is an irritant. Contact with the skin, mucous membranes or eyes should be avoided. If case of skin contact, immediately wash with soap and water. If the material should come into contact with the eyes, immediately rinse with copious amounts of water for at least 15 minutes and consult a physician or ophthalmologist. In rare cases, contact with the skin may lead to a sensitization towards methacrylates. Commercial medical gloves do not provide protection against the sensitizing effects of methacrylate.

Storage Conditions : Always store in dry place between 10°C to 24°C. Close bottle with the original cap immediately after use. Close the Container carefully after use. Keep Away from Moisture. Shelf Life not affected after opening if stores in a cool dry place. If not to be used for a longer duration, keep them in refrigerator. Allow material to reach room temperature prior to use.

Shelf Life : 3 years from the date of manufacturing.

Presentation :

Ref : 10001

Intro Pack : 5ml Plastic Dropper Bottle in a mono carton.

Ref : 10002

Economy Pack : 2 x 5ml Plastic Dropper Bottles + 40 x Micro Brushes in a kit

Disposal conditions : To be disposed off based on the prevailing rules of the country in which it is used.

Prevest DenPro Limited
Unit I, 38, Industrial Estate,
Digiana, Jammu-180010, India
Unit II, EPIP, Bari Brahmna,
Jammu-181133, India
www.prevestdenpro.com
Customer Care No. : +91 1923 222779 (10AM-5PM)
Email : info@prevestdenpro.com

EC REP
Obelis s.a
Boulevard General Wafis 53
1030 Brussels, Belgium
Tel: +(32) 2 732 59 54, Fax: +(32) 2 732 60 03
E-mail : mail@obelis.net

CE
0123

Revision :	02
Date :	27-02-2017

LOT Symbol for "BATCH CODE" **USE BY** Symbol for "USE BY" **NON STERILE** Symbol for "NON STERILE" **KEEP DRY** Symbol for "KEEP DRY"

CAUTION Symbol for "CAUTION" **CONSULT INSTRUCTION FOR USE** Symbol for "CONSULT INSTRUCTION FOR USE" **REF** Symbol for "CATALOGUE NUMBER"

EC REP SYMBOL FOR EUROPEAN COMMUNITY AUTHORIZED REPRESENTATIVE **DO NOT USE IF PACKAGE IS DAMAGE** Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"

KEEP OUT OF THE REACH OF CHILDREN Symbol for "KEEP OUT OF THE REACH OF CHILDREN" **MANUFACTURER NAME & ADDRESS** Symbol for "MANUFACTURER NAME & ADDRESS"

TEMPERATURE LIMITATION Symbol for "TEMPERATURE LIMITATION" **KEEP AWAY FROM SUNLIGHT** Symbol for "KEEP AWAY FROM SUNLIGHT" **FLAMMABLE** Symbol for "FLAMMABLE"

CORROSIVE Symbol for "CORROSIVE"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding.