

MTA Plus



Rating: Not Rated Yet

Price

Variant price modifier:

?. 2495.00

Price with discount ?. 2495.00

Sales price ?. 2495.00

Discount

[Ask a question about this product](#)

Description

Description:

MTA Plus is a specialized dental cement based on Di and Tri-Calcium Silicate, derived from advanced material research in inorganic hydraulic powder technology. MTA Plus is indicated for both vital pulp therapy and Endodontics procedures.

Indications:

- .Vital pulp
- . Cavity Lining
- . Base
- . Pulp-capping
- . Pulpotomies
- Non-Vital / endodontics*
- . Obturation
- . Sealing
- . Apexification
- . Root-end filling
- . Perforation healing
- . Root Resorption

Product Benefits:

- .The original bioactive cement formula
- . Finer powder for smoother, easier mixing, handling and placement
- . Handy packaging- protects powder in desiccant-lined container
- . Dispense powder as needed
- . Optional gel for stability and washout resistance.

Endodontics: MTA Plus

- . Non-cytotoxic, anti- microbial in vitro against common bacteria and E. Faecails
- . No special equipment or dispensing instrument needed

Product Packaging:

Trial Pack: 1 x Jar 1g, 1 x Gel Bottle, 5ml

Video - <https://www.youtube.com/watch?v=9ZCDzki0Y4U>